



IVetMATEIV Fluid Warmer

- Easily hangs from an IV pole or cage
- •Slim design and rubber feet allow unit to lay flat
- Strong magnets hold door closed
- 10 foot power cord

AIV Part #	Part Description
VM13800	IVetMATE Fluid Warmer
AC13761	Replacement Power Cord

Reliable & Accurate Fluid Warming at Multiple Distances & Flow Rates

The IVetMATE is for Veterinary Use Only.

Optimal Distance from Exit of the Warmer (T = 100 - 98°F / 37° - 36°C)		
Inches		
4.5-5		
6.5-8		
10-11		
12-14		
15-17		
17-20		
16-20		

aiv-vet.com • 800.516.1435 sales@aiv-vet.com

